

What is your biosensor good for?



A model to predict the accuracy
of biochemical measurements
made with two-state sensors

Julian Stanley

Apfeld Lab, Northeastern University


International *C. elegans* Meeting, June 2019

We use fluorescent biosensor proteins to make absolute measurements in cells




Absolute measurements allow us to make predictions that relative measurements do not



Relative measurements




Absolute measurements




We use biosensors with two states




We use a mathematical map to make absolute measurements in cells




Errors in R influence redox potential measurements




Errors in R influence redox potential measurements




Errors in R influence redox potential measurements




Errors in R influence redox potential measurements





Errors in R influence redox potential measurements




Errors in R influence redox potential measurements




We can estimate the error associated with a true redox potential




How accurately can we measure E_{GSH} ?




How accurately can we measure E_{GSH} ?




A phase plot shows how accurately we can measure E_{GSH}




A phase plot tells us what E_{GSH} values our sensor is well-suited to measure




A phase plot tells us what E_{GSH} values our sensor is well-suited to measure




A phase plot tells us what E_{GSH} values our sensor is well-suited to measure





A phase plot tells us what E_{GSH} values our sensor is well-suited to measure



E_{GSH} values our sensor is well-suited to measure




Suitable ranges can be compared within sensor families



* Based on empirically-derived spectra characteristics

The Sensor Overlord can tell you what your favorite sensor can measure accurately (<https://sensoroverlord.org>)



S.O. Home Custom Sensor Settings

Upload Spectra Input Characteristics

Upload spectra as .csv

Browse... No file selected

Sensor Type

Redox

pH

Other

Midpoint Value (E0, pKa, etc) Note: Ignored for 'other' sensor type

1

This screenshot shows the "Input Characteristics" tab of the application. It includes fields for uploading spectra (.csv) and selecting a sensor type (Redox, pH, or Other). A midpoint value field contains the number 1. The "Sensor Type" section has a radio button for "Redox" which is currently selected.

Discussion

We can now can:

- Know which values of E_{GSH} we can measure accurately in the pharynx
- Quantify accuracy improvements from better imaging techniques
- Compare the accuracy of two-state sensors
- Identify of what new sensors are needed
- Reclaim underused sensors

Thank You!



Javier
Apfeld
(Fearless
Leader)



Jodie
Norris
(1014B)



Frank
Servello
(1008B)



Sean
Johnsen
(300B)





And a special thanks to Erin Cram's lab, our Northeastern worm buddies!
Avery Lord (240B), Charlotte Kelley (241C), Hannah Pettit (428A), Perla
Castaneda (423B), Jeff Bouffard, and others!



Questions?

Appendix



How accurate is a biosensor?




How do we tell the difference between quantitative sensors?

- Kinetics – well-established
- Suitable ranges—based loosely on midpoint potential, pKa, affinity, etc. Not well characterized. Here we characterize it for most biosensors with two-states, if you're taking ratiometric measurements.


Other groups have developed redox biosensors to make measurements in different environments

| Sensor | Design | Midpoint |
|---|---|--------------------|
|  roGFP1-iE | Designed for “more oxidizing” environments | -236 mV |
|  roGFP1-R12 | Designed for best kinetics in moderate environments | -265 mV |
|  Other roGFPs | Designed for a wide variety of environments | -287 mV to -229 mV |


Challenge: how accurate is each biosensor?




We can present the analysis results as a phase plot



The phase plot depends on microscope precision




By predicting accuracy, we created a framework to compare similar biosensors




Map between spectra and biochemistry depends on the sensor and ratio wavelengths



Specific examples



General case



We observe a ~2.8% relative error in microscopy



The framework generalizes to most two-state biosensors

Biochemical Property

$$E_{GSH} = E'^\circ - \frac{RT}{2F} \ln\left(\frac{\text{Reduced}}{\text{Oxidized}}\right)$$

$$pH = pK_a - \log\left(\frac{[HA]}{[A^-]}\right)$$

Generally:

$$\text{Biochemical Property} = C - D \log_b\left(\frac{\text{State}_2}{\text{State}_1}\right)$$

State Probability

$$\frac{[\text{Oxidized}]}{[\text{Oxidized}] + [\text{Reduced}]} = \frac{R - R_{\text{Red}}}{R - R_{\text{Red}} + \delta_{\lambda_2} (R_{\text{Ox}} - R)}$$

$$\frac{[A^-]}{[A^-] + [HA]} = \frac{R - R_{[HA]}}{R - R_{[HA]} + \delta_{\lambda_2} (R_{[A^-]} - R)}$$

Generally:

$$\frac{[\text{State}_1]}{[\text{State}_1] + [\text{State}_2]} = \frac{R - R_{\text{State}_2}}{R - R_{\text{State}_2} + \delta_{\lambda_2} (R_{\text{State}_1} - R)}$$